

EXTENDED EQUIPMENT SHUTDOWN PROCEDURE

ENSURE EQUIPMENT RELIABILITY WHEN YOU RETURN

The extended shutdown of equipment requires additional effort on top of the regular end-of-day disinfection procedure to avoid maintenance issues.

DAILY & WEEKLY SUCTION MAINTENANCE



SMART Suction Maintenance Videos can be found on our website: cattani.com.au/videos/



WHY IS THIS IMPORTANT?

During extended breaks the turbine in the suction motor can seize, causing maintenance issues and delays when you return.

The suction turbine has very fine clearance tolerances and any humidity left in the system can cause a thin oxidative build-up during an extended break. Many clinics come back from breaks to find their equipment won't start.

With some extra care taken during shutdown this can be avoided.

KEY POINTS OF FINAL SUCTION SHUTDOWN

1. Ensure the suction pipework is clear of fluid prior to disinfection. This will ensure that the disinfectant solution is not diluted in the pipework.
2. Ensure the PULI-JET Gentle 2.0 disinfectant solution is mixed to the correct concentration. 30ml of PULI-JET Gentle 2.0 disinfectant to 2.5 litres of water (per chair).
3. Aspirate the solution using a pulsed action. Pulse cleaners perform this task automatically.
4. After completing the disinfection process, ensure all suction hoses are activated and aspirating dry air for minimum 1 hour to ensure residual moisture at the suction unit is removed.
5. Suction and compressors need to have mains power supplies switched off at the powerpoint. We also recommend you unplug the unit from the wall

RESTARTING EQUIPMENT AFTER AN EXTENDED SHUTDOWN

When planning to reopen for business, we highly recommend that clinics run through tests with all equipment before commencing to see patients.

We recommend doing this a few days before so you will have time to arrange any technical support should you need it.

Don't wait until the day of reopening and find that something is not operating properly.



Australia

+61 3 9484 1120
cattani@cattani.com.au
www.cattani.com.au
280 Dundas Street
Thornbury VIC 3071

New Zealand

0800 68 22 88
cattani@cattani.co.nz
www.cattani.co.nz