

E COMPRESSORS

PLEASE NOTE:

- This is a manual supplement only – **please read the accompanying manual published by Cattani S.p.A. supplied with the unit prior to use** (contact Cattani Australia if you do not have it).
- This unit can only be installed, maintained or worked on by a trained and competent technician.

This compressor is fitted with an SMC condensate trap, which is recommended for humid environments.

Function:

This condensate trap removes the bulk of moisture from compressed air, relieving the load on the main drying system.

Connection to drain:

The drain from the SMC condensate trap (1) must be connected to a waste point within the plant-area, such as a tundish, in accordance with all local plumbing regulations.

The waste point must be lower than the auto-drain on the SMC condensate trap.

Warning:

Failure to connect to a drain point will result in moisture and fluids collecting underneath and around the compressor.

Please phone the Cattani Tech Support Team on 03 9484 1120 if you would like any further information.



1

Doc Name: MSUP001 E Compressors	V. No: 1.1	Expiry date: 17/09/2014	Content Owner: QM
Issue Date: 17/09/2012	Updated by: Andrew	Updated: 3/9/13	Print date: 18/09/2019
Authorised by: QM			

Updates or changes to this document may not be made without the express permission of the Content Owner.