

Blok Jet Silent 2

Intended purpose:

Blok Jet Silent 2 is a semi-wet suction system and dental compressor intended for use in dental clinics.

PLEASE NOTE:

- The Blok Jet Silent 2 must be installed inside the clinic, with an ambient temperature of +5 to +25 degrees celcius.
- This is a manual supplement only – **please read the accompanying manual published by Cattani S.p.A. supplied with the unit prior to use** (contact Cattani Australia if you do not have it).
- This unit can only be installed, maintained or worked on by a trained and competent technician.

Installation – Electrical:

- This unit should be supplied by 2 X dedicated 240V General Power Outlets (GPO), which are protected by a Residual Current Device (RCD).
- This unit should be powered off when the clinic is not in use (i.e. power to the circuit should be isolated off overnight).
- Please note that there are two power supplies – both must be isolated before any maintenance or work is carried out (and by trained and competent persons only).

Installation – Drain Connection:

- The drain of the Turbo SMART must be connected to a drain as per the requirements of AS/NZS 3500.
- The drain the Turbo SMART connects to must be connected to waste via Watermarked Water Seal (i.e. trapped drain that carries the Watermark logo).
- The drain hose, hose clamps, and tube end (included in the 'plumber's kit' supplied with this unit) must be used to connect the fluid outlet to drain. All fittings must be tightened and secured. We recommend that the threaded plastic fittings be tightened with Teflon tape.
- Once the unit has been installed, aspirate no less than 6 litres of water through the surgery suction and observe the unit and connection to drain for any leaks. Any sign of leaks must be remedied so that the connections are sound and have no leaks.

Installation – Exhaust Connection:

- The exhaust connection must be run away from the plant-area, the air-inlet of the air-compressor/s (and any other equipment that has air inlets) and any habitable area, as per "ISO TC 22595-1 2006 Dentistry Plant Area Equipment Part 1 Suction Systems".
- If the above cannot be achieved, a bacterial filter must be fitted to the exhaust air.
- The exhaust air line should be protected from the ingress of fluids (i.e. dog leg on outlets to protect against rain entering).

For information on the Cattani Warranty, please consult the Warranty Statement included on the outer packaging of this unit, or download it from www.cattani.com.au/support.php by selecting 'warranty information' in the 'Choose Category' dropdown menu. Alternatively, you can phone Cattani Australia on 03 9484 1120 for a copy to be sent to you.

Please phone the Cattani Tech Support Team on 03 9484 1120 if you would like any further information.

Doc Name: MSUP0012 Blok Jet Silent	V.No: 1.0	Expiry date: 17/06/2016	Content Owner: QM
Issue Date: 9/7/14	Updated by: Andrew	Last update: 9/7/14	Print date: 9/07/2014
			Authorised by: QM

Updates or changes to this document may not be made without the express permission of the Content Owner.