

# CATTANI



## Suction Tips

### Your In-Surgery Solutions.

Our suction tips are designed to work in harmony with clinic staff and their operating techniques. The range of tips and adaptors give you flexibility to choose the most appropriate tip for each procedure.

The large high volume tips should be used for any procedures that create a spray or aerosol. Germs and bacteria in the spray plume can be inhaled by the dentist and nurse and may lead to possible health issues.

Patient comfort is also very important. Particles, saliva and water can be removed easily from the back of the mouth therefore allowing for uninterrupted treatment sessions.

The small tips are suited to surgical procedures for precise evacuation of fluid from the oral cavity.

To find out which unit suits your surgery's individual needs, head over to [www.cattani.com.au](http://www.cattani.com.au)



# Suction Tips and Adaptors

## Surgical Tips for 11mm hoses

Tips come in packets of three.



**CT9**  
Surgical Tip

Tip width 6mm.  
11mm-d

**CT17**  
Spray Interceptor

Tip width 18mm including flange.  
11mm-d

**CT21**  
Surgical Tip

Tip width 8mm.  
11mm-d

## High Volume Tips for 17mm hoses

Tips come in packets of three.



**CT11**  
High Volume Tip

Tip width 24mm including flange.  
17mm-d

**CT20**  
High Volume Tip

Tip width 15mm including flange.  
17mm-d

**CT22**  
High Volume Tip

Tip width 18mm including flange.  
17mm-d



**Funnel spittoon**

Fits into either the 11mm or 17mm hose terminal.



**Adaptor for saliva ejector**

Adapts most saliva ejector tips to either the 11mm or 17mm hose terminal.



**Reducer**

Reduces the internal bore diameter of the 17mm hose terminal to 11mm for fitment of 11mm suction tips.

All tips and adaptors are autoclavable to 135 degrees.



**Puli-Jet Gentle**

Daily disinfection of your suction system is vital to keep your suction lines clean and maintain your suction unit.

Prevents blockages and unpleasant odours.