

# CATTANI



## Vento

### Limited Space - We Have The Solution...

#### Compact

The Vento System solves the problem of limited space within your plantroom. The frame structure of the Vento System allows the installation of your oil free, dry air compressor above the SMART Suction Unit. This effectively halves the floor space required to house the two units.

#### Benefits

Many combinations are available. If a particular combination is not listed, we can build a special one to suit your requirements upon request. If reduced noise levels are required, we can simply add a 'Q' (Quiet) series compressor to the combination.

#### Why Cattani?

- Minimal installation required
- High vacuum performance
- Consistently deliver clean, dry, oil free hygienic air
- Low maintenance

To find out which unit suits your surgery's individual needs, head over to [www.cattani.com.au](http://www.cattani.com.au)





*Compressors are available with acoustic hoods.*

## Product Range:



**AC200 & Turbo SMART A  
with Amalgam Retention**

Recommended for 2 surgeries

Dimensions: H=1400 W=700 D=540



**AC400 & Tandem Turbo SMART A  
with Amalgam Retention**

Recommended for 5 surgeries

Dimensions: H=1450 W=1155 D=550



**AC300 & Turbo SMART B  
with Amalgam Retention**

Recommended for 3-4 surgeries

Dimensions: H=1460 W=700 D=540



**AC600 & Tandem Turbo SMART B  
with Amalgam Retention**

Recommended for up to 8 surgeries

Dimensions: H=1590 W=1320 D=590