


SAVE YOUR ENERGY

BE SMART ABOUT YOUR SUCTION AND REDUCE YOUR POWER BILLS

Cattani SMART Suction Systems are the brains behind your practice, using intelligent technology and design to provide consistent, high-performance suction while cutting energy costs.





Economical – With SMART suction, energy consumption drops significantly – and so do your energy costs – thanks to Cattani-pioneered inverter technology. There is no wasted energy in the form of extra power consumption and heat, and no unwanted fluctuations in air flow common to other systems.

Flexible – There is a system configuration for every clinical scenario, from a single chair up to hospital and university set-ups. “Plug and play” parallel installation also provides built-in redundancy to protect surgery workflow. With inverters monitoring suction demand in harmony, running multiple parallel systems has never been so simple.

Smart – The inverter control senses exactly the real-time requirements of your practice, from the number of chairs in use, right down to size of suction tip and amount of fluid being aspirated, and adjusts the speed and output of the unit accordingly to ensure a consistent air flow. Pick up the suction tip and you’ll see what we mean.

3-year warranty – we stand by our product, and nothing says reliability like an industry-leading warranty.

Simple, SMART suction for every situation

NEW

SMART Cube – Micro SMART and Turbo SMART with acoustic housing. Cuts noise level in half and delivers further energy savings

- **Micro SMART**
– for 1-2 surgeries*
- **Turbo SMART A**
– for 2-3 surgeries*
- **Turbo SMART B**
– for up to 4 surgeries*
- **Upgrade from Turbo SMART A to Turbo SMART B**
– an incredibly simple, non-disruptive process that doubles suction capacity without increasing system footprint

- **Multi SMART** – parallel systems for up to 20 simultaneous users
- **Maxi SMART** – parallel systems for more than 20 simultaneous users
- **Compact solutions for multiple purposes**
– Vento, Blok-Jet Silent

*Tandem options that provide twice the suction capacity are also available



SMART Suction from Cattani

- Inverter technology slashes energy consumption
- Intelligent, consistent high performance
- Easy to add parallel units
- Compact, all-in-one design
- Remote alarm capabilities
- Unrivalled technical service back-up
- Generous industry-leading parts and labour warranty
- Easy installation and low maintenance
- Amalgam retention optional
- Cattani After Care Programme



WE LOVE WHAT WE DO.

Have confidence in Cattani

For more than 50 years, Cattani's global business has been dental air technology – no distractions. We stand for quality, performance and reliability with the aim of keeping the dental practice running day after day, year after year.



Choice is the experience at Cattani

Australia

+61 3 9484 1120

cattani@cattani.com.au

www.cattani.com.au

280 Dundas Street
Thornbury VIC 3071

New Zealand

Toll Free 0800 68 22 88

cattani@cattani.co.nz

www.cattani.co.nz