

WORKING AS ONE



MAGNOLIA
CHEMICAL PROFESSIONALS



CATTANI
AIR TECHNOLOGY
— AUSTRALASIA

MAGNOLIA **by CATTANI –** **FOR CONSTANT** **DISINFECTION**

Designed specifically for cleaning, disinfecting and sanitising dental suction systems and surgeries, Magnolia products can be trusted to keep your practice hygienic.

WE LOVE WHAT WE DO.

Hygienic, Safe, Reliable

Led by Puli-Jet Gentle, Magnolia cleans and protects your suction system from the inside out, for everyone's sake. Make it your suction and surgery disinfecting range of choice.



 **MAGNOLIA**
CHEMICAL PROFESSIONALS



Suction maintenance at your fingertips – just point your smartphone camera at this QR code to view our short instructional videos, or go to www.cattani.com.au/videos



Puli-Jet Gentle and Pulse Cleaner

Daily suction system maintenance starts with Puli-Jet Gentle and the Pulse Cleaner agitator. Disinfection takes just a few minutes.

- Ideal for all suction systems
- Commercial grade disinfectant
- Dissolves blood and mucus
- Removes scaling from pipes and tubing
- Leaves behind an anti-scaling surfactant
- 1-litre bottle enough for 33 cleaning cycles
- Pulse Cleaner automatically creates turbulence for optimal effectiveness



Antifoaming Disinfectant Tablets

Foam created by blood and mucus can cause problems in suction systems. An Antifoaming Disinfectant Tablet placed in the filter basket each day dramatically reduces this foam.

- No more foam issues
- Slow disinfectant release for long-lasting effectiveness
- Prevents need for harsher, more damaging cleaning products
- Pack contains 50 tablets; one per chair per day



Lubri-Jet Silicone Drops

Filter o-rings need lubricating, and Lubri-Jet is ideal for the task.

- Preserves o-ring seal and prolongs ring life
- Use daily on chairside filter o-ring
- Use weekly on suction unit filter o-ring
- Useful on instruments with friction-prone parts before sterilisation
- 4 x 30ml bottles per pack



Eco-Jet 1

Disinfect your surgery surfaces as required with this commercial grade spray.

- Low alcohol content, pleasant odour
- Comes in ready-to-use 500ml bottle

Magnolia Works For You

Correct and regular cleaning with Magnolia protects the system and reduces deterioration that can impede suction performance.



Harmony happens at Cattani

Australia

+61 3 9484 1120

cattani@cattani.com.au

www.cattani.com.au

280 Dundas Street
Thornbury VIC 3071

New Zealand

0800 68 22 88

cattani@cattani.co.nz

www.cattani.co.nz