



SAVE THAT IMPORTANT FILTER CHANGE DATE

Dear Doctor

Thank you for choosing a Cattani compressor.

Behind the design of this machine is more than 40 years of research and development as well as Cattani's unshakeable commitment to providing you with clean, dry air. This compressor will reward you with a long and reliable life if placed in the correct position, with a good supply of clean air and with minimal maintenance as well as periodic checks by a qualified technician, using original spare parts.

The compressor head has been fitted with standard filters. In this pack you will find **an adhesive label on which to write the filter product code and operation start date**. The label should be attached to the compressor in a clearly visible position. In 12 months when you order new filters, the replacement package will contain a new label, which you can place over the first one. **The item code for a standard filter is K-168563. The code for an anti-bacterial HEPA filter is K-168429.**

Please note:

With clean air, the filters will last 12 months. With dust-rich air, check the filters frequently and act accordingly.

Should you need further help or information, don't hesitate to contact our service team.

Yours sincerely

Cattani S.p.A.