



# YOUR QUIET SPACE

Put your space and noise issues behind you with BLOK-JET SILENT 2 – your amazingly compact, quiet plant room solution.

Takes up less than 0.5 m<sup>2</sup>

Serves up to four chairs

Turbo SMART B suction unit  
+ AC300 compressor



**Blok-Jet Silent 2 is your plant room.  
Contact your Cattani dealer now.**



Find your quiet space with Cattani  
[www.cattani.com.au](http://www.cattani.com.au)

## Blok-Jet Silent 2 – the self-contained plant solution for up to four chairs

The Blok-Jet Silent 2 solution is to stack the units inside their own cabinet. It can service up to four chairs and is ultra-quiet, so is an excellent option when mitigating noise is a priority. The cabinet has a very small footprint and stands on caster wheels, so is easy to service.

- Compact – the size of a fridge
- Ultra-quiet
- Mobile – it stands on lockable wheels
- Cabinet doors also lockable
- External display panel shows compressor and suction status

### BLOK-JET SILENT 2

<b>NO. OF CHAIRS</b>	4
<b>HEIGHT MM</b>	1935
<b>WIDTH MM</b>	760
<b>DEPTH MM</b>	650
<b>WEIGHT KG</b>	270
<b>SOUND OUTPUT dB(A)</b>	54



### UNIT INCLUDES

AC300 & Turbo SMART B with ISO 18 Hydrocyclone



### AC300

<b>NO. OF CHAIRS</b>	3-4
<b>POWER</b>	1ph 230V 50Hz 10.2A 1.5kW 1ph 220V 60Hz 12.2A 1.75kW
<b>HEIGHT MM</b>	750
<b>WIDTH MM</b>	620
<b>DEPTH MM</b>	520
<b>WEIGHT KG</b>	64
<b>SOUND OUTPUT dB(A)</b>	73.6
<b>MAX AIR FLOW AT 5 BAR NL/min</b>	238
<b>TANK SIZE (LITRES)</b>	45

\* 3ph available  
\* w/out dryer available



### TURBO SMART B

<b>NO. OF CHAIRS</b>	4
<b>MAX AIR FLOW l/min</b>	1750
<b>MAX. VACUUM FOR CONTINUOUS OPERATION mBar</b>	210
<b>POWER</b>	1ph 9A 1.5kW
<b>HEIGHT MM</b>	630
<b>WIDTH MM</b>	390
<b>DEPTH MM</b>	350
<b>WEIGHT KG</b>	39.5
<b>SOUND OUTPUT dB(A)</b>	73.7

**ISO 18 HYDROCYCLONE  
AMALGAM RETENTION**

Optional