


THE PERFECT FIT

Halve your plant
equipment footprint
with **VENTO** –
your vertical plant
room solution.



- Compact
- Flexible design
- Quiet options
- Serves up to 8 chairs*

Vento can incorporate:

- Micro SMART or Turbo SMART suction unit
- Or a quiet SMART Cube suction unit
- Any Cattani compressor, including CAD/CAM
- Quiet compressor hood
- Amalgam retention

* We can increase the number of chairs to suit you – talk to us or talk to your dealer.



Vento – the ultimate space saver for up to eight chairs

Designed and developed by Cattani Australia, the Vento approach is to stack the suction unit and the compressor on a special frame, halving the floor space. This is the perfect plant configuration for practices in shopping centres where space is tight and a dedicated plant room is a luxury. Instead, the Vento can easily fit into a ventilated cupboard.


AC200Q

NO. OF CHAIRS	2-3
SOUND OUTPUT dB(A)	63

MICRO SMART CUBE

NO. OF CHAIRS	1-2
SOUND OUTPUT dB(A)	59

HEIGHT MM	1680
WIDTH MM	695
DEPTH MM	600



ALTERNATIVE VENTO OPTIONS



	SMART Vento with C1A & AC200	SMART Vento Trio with C2A & AC300	SMART Vento Duo Tandem with C1A & AC400	SMART Vento Trio Tandem with C2A & AC600
	DUO	TRIO	DUO TANDEM	TRIO TANDEM
ITEM CODE	C-VC1A20	C-VC2A30	C-VTC1A4	C-VTC2A5
COMPRESSOR	AC200	AC300	AC400	AC600
SUCTION	TURBO SMART A	TURBO SMART B	TANDEM TURBO SMART A	TANDEM TURBO SMART B
NO. OF CHAIRS	2	3-4	5	UP TO 8