

# MULTI-SMART PLANT ROOM SOLUTIONS

**NEW** COMPACT DESIGN WITH NO COMPROMISE ON PERFORMANCE





### SMART THINKING, THE WAY WE INNOVATE.

The NEW Multi-SMART suction systems are specifically designed for the Australasian region. Co-engineered, our Italian and Australian engineering and technical teams worked closely to develop a new 12 – 20 chair solution that is more compact with no compromise on performance.

#### Compact design

 Space is a premium, the NEW Multi-SMART has been designed to be more compact than its predecessor.

## Consistent and reliable air flow

 Multi-SMARTS are assembled in parallel, and thanks to inverter technology provide steady and reliable air flow, with built-in redundancy. Our modular design allows for easy HEPA filter upgrade onsite.

# Easy maintenance, minimal disruptions

- Installation of Multi-SMARTS, whether a T3 (12 chairs), T4 (16 chairs) or T5 (20 chairs), is a simple process, the system is designed with easy access for future annual servicing.
- Each suction module of the system can be isolated to allow for easy servicing while still allowing uninterrupted suction, keeping the dental site running without disruption.













MULTI-SMART T3		<b>MULTI-SMART T4</b>	MULTI-SMART T5
ITEM CODE			
Multi-SMART	035105BX3	035105BX4	035105BX5
Multi-SMART with H14 Filtration	035105BHX3	035105BHX4	035105BHX5
Multi-SMART with Hydrocyclone Amalgam Retention	035115BX3	035115BX4	035115BX5
Multi-SMART with Hydrocyclone Amalgam Retention and H14 Filtration	035115BHX3	035115BHX4	035115BHX5
NO. OF CHAIRS	12	16	20
POWER 4 x 1PH 2	40V 9A 1.87 KW	5 x 1PH 240V 9A 1.87 KW	6 x 1PH 240V 9A 1.87 KW
HEIGHT MM	1280	1280	1280
<b>WIDTH</b> MM 1000 (1030 with	h control panel)	1345 (1375 with control panel)	1690 (1720 with control panel)
DEPTH MM	760	760	760
WEIGHT KG	148	184	220
SOUND OUTPUT dB(A)	78.5	79.7	80.7
MAX AIR FLOW L/min	5100	6800	8500
MAX AIR FLOW FOR CONT. OPERARTION mBar	210	210	210

# Sales and Technical Support:

Australia +61 3 9484 1120 cattani@cattani.com.au www.cattani.com.au New Zealand 0800 68 22 88 cattani@cattani.co.nz www.cattani.co.nz

### **Company Address**

280 Dundas Street, Thornbury VIC 3071

