CATTANI COMPRESSORS

AC100 1 CHAIR

- Clean, Dry, Oil-Free Air
- 7-year Parts and Labour Warranty
- Outstanding Performance and Reliability
- Sound Reduction Options

Reliable: Robust Cattani compressors keep practices running, delivering clean, dry, oil-free air with such a level of reliability that we give it a 7-year warranty. The AC100, for 1 chair, is easy to install and low maintenance.

Quality Protection: The tanks in all Cattani compressors have a patented silver ion coating that delivers highly effective protection against microbes such as bacteria and mould. Permanent antimicrobial protection is assured. Internal components are also specially coated to reduce mechanical wear and tear.

Sound Reduction: If the compressor is in close proximity to treatment and staff rooms, Cattani offers the AC100Q, which significantly reduces noise output.





PERFORMANCE AND RELIABILITY

Cattani Australasia is a family company that stands for outstanding quality, performance and service. We have more than 30 years' experience in the field of dental air technology, and are proud to work with what we believe are the best compressors and suction systems in the business.

The confidence we have in our compressors' superior design, mechanical precision and robust components means we have no hesitation in backing them with an industry-leading 7-year parts and labour warranty.

Cattani compressors, which come in sizes to suit every clinic, are manufactured to the strictest standards in Parma, Italy. Each unit is rigorously tested before delivery, allowing the clinician to have the utmost confidence in its performance and reliability.

SOUND REDUCTION OPTIONS

In some clinics, plant room equipment noise can be an issue, for example where the plant room is close to treatment and staff rooms. In these situations, we have sound reduction options that can dramatically cut noise.









AC100

NO. OF CHAIRS	1
POWER	1~ 230V 50Hz 3.8A 0.55kW
HEIGHT MM	720
WIDTH MM	600
DEPTH MM	470
WEIGHT KG	47
SOUND OUTPUT dB(A)	70
MAX AIR FLOW AT 5 BAR NL/min	67.5
TANK SIZE (LITRES)	30
PRODUCT CODE	(AU/NZ) K-AC1000, (IT) 13125











AC100G

NO. OF CHAIRS	1
POWER	1~ 230V 50Hz 3.8A 0.55kW
HEIGHT MM	860
WIDTH MM	600
DEPTH MM	510
WEIGHT KG	57
SOUND OUTPUT dB(A)	63
MAX AIR FLOW AT 5 BAR NL/min	67.5
TANK SIZE (LITRES)	30
PRODUCT CODE	(AU/NZ) K-AC100Q, (IT) 13144



ACOUSTIC COVER

H = 887 mm W = 870 mm D = 700 mm

Reduces noise of AC100 compressor to 50.7 dB(A)

Product code: K-010800



VENTILATION KITS

Recommendations – 100 m3/h for each kW of equipment installed. 250 & 500 kits include fan, ducting and vent cover. 800 includes fan only.

Product codes:

250 – L-TD2500-KIT 500 – L-TD5000-150

800 - L-TD8000



HEPA H14 FILTER

For medical quality air. Non-autoclavable.

Product code: K-168429



DRYER PRE FILTER FOR COMPRESSOR (FOR HUMID CLIMATES)

A dryer pre filter is recommended for humid climates. This additional moisture trap is fitted prior to the drying column. The pre filter has proven extremely effective and reliable in increasing the efficiency of the drying system by removing bulk condensation from the air prior to entering the drying column. Can be retro-fitted.

Product code: K-E10000

Australia

+61 3 9484 1120 cattani@cattani.com.au www.cattani.com.au

280 Dundas Street Thornbury VIC 3071

New Zealand

0800 68 22 88 cattani@cattani.co.nz www.cattani.co.nz

