

AFTER CARE TRAINING PROGRAMME™ ANNUAL SERVICING SCHEDULE

AC & CADCAM COMPRESSORS

AFTER CARE
∴ 100% RELIABILITY



AC RANGE	AC100	AC100Q	AC200	AC200Q	AC300	AC300Q	AC4001P	AC400Q1P	AC400	AC400Q	AC600	AC600Q	AC900	AC1200	AC1800
CADCAM			AC210	AC210Q	AC310	AC310Q	AC4101P	AC410Q1P	AC410	AC410Q	AC610	AC610Q	AC910		
YEAR															
1	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
2	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
3	1 Kit K-SK101B	1 Kit K-SK101B	1 Kit K-SK201B	1 Kit K-SK201B	1 Kit K-SK301B	1 Kit K-SK301B	2 Kits K-SK201B	2 Kits K-SK201B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
4	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
5	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
6	1 Kit K-SK101B	1 Kit K-SK101B	1 Kit K-SK201B	1 Kit K-SK201B	1 Kit K-SK301B	1 Kit K-SK301B	2 Kits K-SK201B	2 Kits K-SK201B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
7	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
8	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
9	1 Kit K-SK101B	1 Kit K-SK101B	1 Kit K-SK201B	1 Kit K-SK201B	1 Kit K-SK301B	1 Kit K-SK301B	2 Kits K-SK201B	2 Kits K-SK201B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B
10	1 Kit K-SK1013B	1 Kit K-SK1013B	1 Kit K-SK2013B	1 Kit K-SK2013B	1 Kit K-SK3013B	1 Kit K-SK3013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK2013B	2 Kits K-SK3013B	2 Kits K-SK3013B	3 Kits K-SK3013B	2 Kits K-SK0106-11B	3 Kits K-SK0106-11B

PLEASE NOTE: ALL SERVICING SCHEDULES ALSO APPLY TO THE COMPRESSOR WITH PRE-FILTER (AC100 "E" VERSION E.G. AC100E).

MINOR SERVICE KITS (Filters are replaced annually to maintain performance)

Item No.	Includes:
K-SK1013B	1x Air Inlet Filter for Compressor Head Unit, 1x HEPA H14 Drying Column Filter & 1x Drying Column O-ring Seal
K-SK2013B	2x Air Inlet Filter for Compressor Head Unit, 1x HEPA H14 Drying Column Filter & 1x Drying Column O-ring Seal
K-SK3013B	3x Air Inlet Filter for Compressor Head Unit, 1x HEPA H14 Drying Column Filter & 1x Drying Column O-ring Seal
K-SK0106-11B	6x Air Inlet Filter for Compressor Head Unit, 2x HEPA H14 Drying Column Filter & 2x Drying Column O-ring Seal

INTERMEDIATE SERVICE KITS (1 PH Compressor Only - Performed Every 3 Years)

Item No.	Includes:
K-SK101B	1x Air Inlet Filter for Compressor Head Unit, 1x HEPA H14 Drying Column Filter, 1x Drying Column O-ring Seal, 1x 20 µF Motor Start Capacitor & 2x Cable Tie
K-SK201B	2x Air Inlet Filter for Compressor Head Unit, 1x HEPA H14 Drying Column Filter, 1x Drying Column O-ring Seal, 1x 40 µF Motor Start Capacitor & 2x Cable Tie
K-SK301B	3x Air Inlet Filter for Compressor Head Unit, 1x HEPA H14 Drying Column Filter, 1x Drying Column O-ring Seal, 1x 55 µF Motor Start Capacitor & 2x Cable Tie

MINOR & INTERMEDIATE SERVICE KIT PARTS

- 1 Air Inlet Filter for Compressor Head Unit
- 2 HEPA H14 Filter for Drying Column & O-ring Seal
- 3 Motor Start Capacitor & Cable Ties (1 PH Only)



DOWNLOAD NOW
AFTER CARE COMPRESSOR SERVICE REPORT

Point your smartphone camera at the QR code to view.

