

# AFTER CARE TRAINING PROGRAMME™ ANNUAL SERVICING SCHEDULE

## SMART SUCTION

**AFTER CARE**  
∴ 100% RELIABILITY



SMART RANGE	Micro SMART	Micro SMART +Hydro.	Micro SMART CUBE	Micro SMART CUBE +Hydro.	Turbo SMART A & B	Turbo SMART A & B +Hydro.	Turbo SMART A & B CUBE	Turbo SMART A & B CUBE +Hydro.	Tandem Micro SMART	Tandem Micro SMART +Hydro.	Tandem Turbo SMART A & B	Tandem Turbo SMART A & B +Hydro.	Multi SMART T3	Multi SMART T3 +Hydro.	Multi SMART T4	Multi SMART T4 +Hydro.	Multi SMART T5	Multi SMART T5 +Hydro.	
YEAR																			
<b>1</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	2x C-211956 1x M-040727	2x C-211956 1x M-040727	2x C-211974 1x M-040727	2x C-211974 1x M-040727	3x C-211974 1x M-040727	3x C-211974 1x M-040727	4x C-211974 1x M-040727	4x C-211974 1x M-040727	5x C-211974 1x M-040727	5x C-211974 1x M-040727	
<b>2</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	2x C-211956 1x M-040727	2x C-211956 1x M-040727 2x C-MR60SK	2x C-211974 1x M-040727	2x C-211974 1x M-040727 2x C-1H00SK	3x C-211974 1x M-040727	3x C-211974 1x M-040727 3x C-1H00SK	4x C-211974 1x M-040727	4x C-211974 1x M-040727 4x C-1H00SK	5x C-211974 1x M-040727	5x C-211974 1x M-040727 5x C-1H00SK	
<b>3</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	2x C-211956 1x M-040727	2x C-211956 1x M-040727	2x C-211974 1x M-040727	2x C-211974 1x M-040727	3x C-211974 1x M-040727	3x C-211974 1x M-040727	4x C-211974 1x M-040727	4x C-211974 1x M-040727	5x C-211974 1x M-040727	5x C-211974 1x M-040727	
<b>4</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	2x C-211956 1x M-040727	2x C-211956 1x M-040727 2x C-MR60SK	2x C-211974 1x M-040727	2x C-211974 1x M-040727 2x C-1H00SK	3x C-211974 1x M-040727	3x C-211974 1x M-040727 3x C-1H00SK	4x C-211974 1x M-040727	4x C-211974 1x M-040727 4x C-1H00SK	5x C-211974 1x M-040727	5x C-211974 1x M-040727 5x C-1H00SK	
<b>5</b>	1x C-211956 1x M-040727 1x C-MSMSK0	1x C-211956 1x M-040727 1x C-MSMSK0	1x C-211956 1x M-040727 1x C-MSCMSK	1x C-211956 1x M-040727 1x C-MSCMSK	1x C-211974 1x M-040727 1x C-TSMSK0	1x C-211974 1x M-040727 1x C-TSMSK0	1x C-211974 1x M-040727 1x C-TSCMSK	1x C-211974 1x M-040727 1x C-TSCMSK	2x C-211956 1x M-040727 2x C-MSMSK0	2x C-211956 1x M-040727 2x C-MSMSK0	2x C-211974 1x M-040727 2x C-TSMSK0	2x C-211974 1x M-040727 2x C-TSMSK0	3x C-211974 1x M-040727 3x C-TSMSK0	3x C-211974 1x M-040727 3x C-TSMSK0	4x C-211974 1x M-040727 4x C-TSMSK0	4x C-211974 1x M-040727 4x C-TSMSK0	5x C-211974 1x M-040727 5x C-TSMSK0	5x C-211974 1x M-040727 5x C-TSMSK0	
<b>6</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	2x C-211956 1x M-040727	2x C-211956 1x M-040727 2x C-MR60SK	2x C-211974 1x M-040727	2x C-211974 1x M-040727 2x C-1H00SK	3x C-211974 1x M-040727	3x C-211974 1x M-040727 3x C-1H00SK	4x C-211974 1x M-040727	4x C-211974 1x M-040727 4x C-1H00SK	5x C-211974 1x M-040727	5x C-211974 1x M-040727 5x C-1H00SK	
<b>7</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	2x C-211956 1x M-040727	2x C-211956 1x M-040727	2x C-211974 1x M-040727	2x C-211974 1x M-040727	3x C-211974 1x M-040727	3x C-211974 1x M-040727	4x C-211974 1x M-040727	4x C-211974 1x M-040727	5x C-211974 1x M-040727	5x C-211974 1x M-040727	
<b>8</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	2x C-211956 1x M-040727	2x C-211956 1x M-040727 2x C-MR60SK	2x C-211974 1x M-040727	2x C-211974 1x M-040727 2x C-1H00SK	3x C-211974 1x M-040727	3x C-211974 1x M-040727 3x C-1H00SK	4x C-211974 1x M-040727	4x C-211974 1x M-040727 4x C-1H00SK	5x C-211974 1x M-040727	5x C-211974 1x M-040727 5x C-1H00SK	
<b>9</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211956 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	1x C-211974 1x M-040727	2x C-211956 1x M-040727	2x C-211956 1x M-040727	2x C-211974 1x M-040727	2x C-211974 1x M-040727	3x C-211974 1x M-040727	3x C-211974 1x M-040727	4x C-211974 1x M-040727	4x C-211974 1x M-040727	5x C-211974 1x M-040727	5x C-211974 1x M-040727	
<b>10</b>	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211956 1x M-040727	1x C-211956 1x M-040727 1x C-MR60SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	1x C-211974 1x M-040727	1x C-211974 1x M-040727 1x C-1H00SK	2x C-211956 1x M-040727	2x C-211956 1x M-040727 2x C-MR60SK	2x C-211974 1x M-040727	2x C-211974 1x M-040727 2x C-1H00SK	3x C-211974 1x M-040727	3x C-211974 1x M-040727 3x C-1H00SK	4x C-211974 1x M-040727	4x C-211974 1x M-040727 4x C-1H00SK	5x C-211974 1x M-040727	5x C-211974 1x M-040727 5x C-1H00SK	

PLEASE NOTE: MULTI SMART RANGE - HEPA FILTERS ARE TO BE REPLACED ANNUALLY IF FITTED (MULTI SMART T3 = 3x C-160540 | MULTI SMART T4 = 4x C-160540 | MULTI SMART T5 = 5x C-160540).

### MINOR SERVICE PARTS

Item No.	Description:
C-211956	1x Micro SMART Suction Inlet Filter O-Ring
C-211974	1x Turbo SMART Suction Inlet Filter O-Ring
M-040727	4x Magnolia Lubri-Jet 30ml Bottle
C-160540	1x HEPA H14 Filter for 4 to 6 Surgeries

### INTERMEDIATE SERVICE KITS (Hydrocyclone Amalgam Retention Systems)

Item No.	Description:
C-MR60SK	ISO 6 Hydrocyclone Service Kit
C-1H00SK	ISO 18 Hydrocyclone Service Kit
<b>Optional:</b>	
C-1H02SK	ISO 18 Hydrocyclone Service Kit (Valve & Kit)

In hard water areas the C-1H02SK can be used instead of the C-1H00SK for Turbo SMART & Turbo SMART CUBE.

### MAJOR SERVICE KITS\*

Item No.	Description:
C-MSMSK0	Micro SMART Major Service Kit
C-MSCMSK	Micro SMART CUBE Major Service Kit
C-TSMSK0	Turbo SMART Major Service Kit
C-TSCMSK	Turbo SMART CUBE Major Service Kit

\*Major service kits contain separate kits for the centrifugal separator, suction motor and internal hoses.

### MINOR, INTERMEDIATE & MAJOR SERVICING

- 1 Suction Inlet Filter O-Ring
- 2 Hydrocyclone Service Kit
- 3 Major Service Kit



**DOWNLOAD NOW**  
AFTER CARE SMART SUCTION  
SERVICE REPORT

Point your smartphone camera at the QR code to view.

